

1. UCSF Chimera, which allows the phi and psi of each residue to be specified.
 Tools -> Structure Editing -> Build -> Peptide
 - paste the sequence in the box
 - assign fully extended phi, psi values to all non-P residues (-179,179)
 - give all Pro residues PPII dihedral angles (-80, 160)
 - save as CTD_chimera.pdb
1. Check the sequence is correct in the PDB file. Check all dihedrals are trans (they are)
 - a. `g_chi -s CTD_chimera.pdb -f CTD_chimera.pdb -omega`
2. `pdb2gmx -f CTD_chimera.pdb -o init.gro -p init.top -ignh -ter -vsite hydrogens -v #CHARMM22STAR`
3. `editconf -f init.gro -o editconf.gro -bt dodecahedron -d 1.5 -center 0 0 0`
4. `genbox -cp editconf.gro -cs spc216 -o genbox.gro -p init.top`
5. `/usr/local/gromacs/4.6/465-mpi403-fftw332-gcc447-2/bin/grompp -f em.mdp -c genbox.gro -p init.top -o genbox.tpr` [make the tpr for genion]
6. `echo 13 | genion -s genbox.tpr -conc 0.15 -o genbox_genion.pdb -neutral -p init.top`

Energy Minimization

7. `/usr/local/gromacs/4.6/465-mpi403-fftw332-gcc447-2/bin/grompp -f em.mdp -c genbox_genion.pdb -p init.top -o em.tpr`
8. `mdrun -s em.tpr`